

# Assessment of drug-induced structural toxicity in human iPSC-derived cardiomyocytes, an image-based approach to complement traditional in vitro safety studies

SP092

Brigitta Rita Szabó<sup>1,2</sup>, Georgios Kosmidis<sup>1</sup>, Paul Volders<sup>2</sup>, Elena Matsa<sup>1</sup>

<sup>1</sup>Ncardia Services B.V., Discovery Technology, Leiden, Netherlands (The), <sup>2</sup>Cardiovascular Research Institute Maastricht (CARIM), Maastricht, Netherlands (The)


## BACKGROUND AND PURPOSE OF THE STUDY

Human induced pluripotent stem cell-derived cardiomyocytes (hiPSC-CMs) are novel alternatives to classic models in the field of safety pharmacology. Significant progress has already been made towards validating the use of hiPSC-CM electrophysiology assays to predict pro-arrhythmic potential. Yet, for other clinically relevant parameters such as structural toxicity – which could facilitate more accurate in vitro safety assessments – further optimization is still warranted.

This study aimed to optimize a cost-effective, scalable method to detect compound-induced morphological changes using high content imaging in combination with MEA recordings and RNAseq to facilitate more accurate *in vitro* safety assessment.

## 1 METHODS

### Experiment design



Multiple lines of hiPSC-CMs (e.g., Ncyte CMs, NCRM-5), manufactured using Ncardia's proprietary protocols, were cultured according to manufacturer's instructions in chemically defined, serum-free medium. hiPSC-CMs were treated for a period of 24 hours with 7 concentrations of each drug, in a range comparable to the clinical maximum plasma concentrations ( $C_{max}$ ) followed by functional, mechanistic and morphological assays.

### Parameters assessed

- Morphology of subcellular target structures
- Functional changes via multi-electrode array (MEA) field potential recordings
- Gene expression analysis (RNAseq)

### Methods

- High-magnification (40x) fluorescence image stacks were acquired using an ImageXpress Micro Confocal microscope. Features such as fluorescence intensity, detected signal area and number of detected objects per cell, were extracted via analysis of 2D image projections using MetaXpress software v6.6
- MEA recordings of electrophysiological alterations including changes in field potential duration (FPD), beat rate (BR) and beat rate irregularity (BRI)
- Mechanism of action of each drug was also investigated via gene expression analysis (RNAseq)

Expected toxicity	Compound	Main mechanism of action
High (A)	Doxorubicin	DNA Intercalator
	Cisplatin	DNA synthesis inhibitor
	Ponatinib	Tyrosine Kinase inhibitor
	Dasatinib	Tyrosine Kinase inhibitor
	Lapatinib	Tyrosine Kinase inhibitor
Intermediate (B)	5-Fluoro-uracil	Antimetabolite
	Methotrexate	DNA synthesis inhibitor
	Omecamtiv Mecarbil	Cardiac myosin activator
	Propofol	GABA <sub>A</sub> receptor modulator
	Bupivacaine	Sodium channel inhibitor
Low (C)	Amiodarone	Class III antiarrhythmic
	Dofetilide	Class III antiarrhythmic
	Digoxin	Na/K ATPase inhibitor
	Chlorpromazine	D2 (dopamine) receptor antagonist
	Erlotinib	Tyrosine Kinase inhibitor
ASA	COX1/2 inhibitor	
Empagliflozin	SGLT2 inhibitor	

Table 1. Compound library

Target structures	Fluorescent probes
Mitochondria	MitoTracker CMX Ros
DNA damage	Anti-γH2A.X antibody
Golgi	Wheat Germ Agglutinin
Endoplasmic reticulum (ER)	Concanavalin A
Gap junctions (connexin 43)	Anti-Cx43 antibody
Lysosomes	LysoTracker Red
Peroxisomes	Anti-PMP70 antibody
Nucleoli	Anti-fibrillarin antibody
Nuclei	DAPI

Table 2. Validated immunofluorescence assays






Image analysis. Fluorescent image (right) and corresponding analysis mask (left) for feature extraction.


## 2 A NON-TOXIC COMPOUND – EMPAGLIFLOZIN



## 3 AMIODARONE – A POTENT ANTIARRHYTHMIC WITH SAFETY RISKS



## 4 HIGHLY TOXIC ANTHRACYCLINE - DOXORUBICIN



## CONCLUSIONS

- Human iPSC-CM derived models are scalable and cost-effective tools for safety assessments
- We were able to optimize staining protocols for a selection of dyes targeting subcellular organelles and membranes, utilizing high content imaging methods, which complement already standardized methods for safety assessment in vitro using hiPSC-CMs
- Using multiple hiPSC-CM lines, our results show the potential of high-throughput image-based structural analysis as a tool in gaining a deeper understanding of morphological changes as a mechanism for drug-induced cardiotoxicity

This project has received funding from the European Union's Horizon 2020 research and innovation program under grant agreement No 858070

