







# Awake & unrestrained PET neuro-imaging

## MADA<sup>©</sup> PET

Motion ADAptive PET





### What is MADA<sup>©</sup> PET?

### Motion adaptive PET imaging solution

- Allows awake and unrestrained brain PET based on motion tracking and motion correction.
- Integrated and easy-to-use methodology
- Minimally invasive head tracking
- Software solution no additional hardware
- Applicable for many types of μPET scanners

# Why MADA<sup>©</sup> PET?


Imaging typically requires anesthesia to avoid stress and/or motion artefacts during imaging.

However, anesthesia may introduce confounding effects such as:


- influence on the normal brain physiology (e.g. perfusion)
- interaction of anesthetic with target molecules
- interaction of anesthetic with radiotracer (e.g. metabolism)
- possible interference with disease development

## How does MADA<sup>©</sup> PET work?


Point source based head tracking using PET data




Motion adaptive reconstruction software



# How does it work in practice?



# Phantom proof of principle



~ 22 MBq [18F]FDG, Siemens Inveon, 10 min scan with motion followed by 20 min static reference scan . Random motion of the phantom was introduced by attaching it to a manually controlled pole.

MADA<sup>©</sup> PET software compensates randomly-introduced motion of the phantom and generates very sharp images with about 1.25mm resolution.

## In vivo proof of principle


Uncorrected image

MADA PET image

MADA<sup>©</sup> software enables awake PET brain imaging

Static reference image

MR reference



~ 37 MBq [18F]FDG, Siemens Inveon, 20 min awake scan followed by 20 min static reference scan

## Unique strengths of MADA PET

- ✓ Enables neuro-imaging without anesthesia
- ✓ Easy-to-use and integrated methodology
  - No additional equipment
  - No additional calibrations
- ✓ Minimally invasive
  - No large and heavy markers
- ✓ Excellent image quality

#### Want to collaborate?



Jeroen Verhaeghe, PhD
Associate professor
jeroen.verhaeghe@uantwerpen.be

Molecular Imaging
Center Antwerp (MICA)
University of Antwerp

MICA: www.uantwerpen.be/en/rg/mica/



Pieter-Jan Guns, PhD
Research and Innovation Manager EGAMI
pieter-jan.guns@uantwerpen.be



EGAMI: www.egami.be